



#### MEMBERS OF THE GROUP

lach:ner

LACHNER, Czech Republic  
[www.lach-ner.cz](http://www.lach-ner.cz)

PROANALYTICA

PROANALYTICA, doo  
[www.proanalytica.com](http://www.proanalytica.com)



REANALLABOR, Hungary  
[www.reanallabor.hu](http://www.reanallabor.hu)

lach:ner

LACHNER Zagreb, Croatia  
[www.lach-ner.hr](http://www.lach-ner.hr)



CENTRALCHEM, Slovakia  
[www.centralchem.sk](http://www.centralchem.sk)

## **LABORATORY FURNITURE, FUME HOODS, VENTILATION SYSTEMS MADE by LACHNER GROUP for YOUR SAFETY and COMFORT**

Following your needs.

Giving you project proposal, design, technically and ergonomically sustainable solutions for your working and storage space.


Offering you comprehensive options for microbiological, chemical, mechanical laboratories and manufacturing facilities.

Delivering top quality, highly chemical resistant and microbiologically acceptable materials.

Installing ventilation, respecting air flow, room/chamber pressure, and gas, liquid and powder exhaust capacity.

Manufacturing and installing according to EU standards.

Type-testing and site testing procedured by accredited laboratories.



## Safety cabinets for safe storage of aggressive, toxic, flammable, water-pollutant substances



90 minutes-fire resistant, highly-acid-proof and corrosion-proof resistant cabinet for safe storage of flammable and corrosive substances



90-minutes-fire resistant cabinet for safe storage of flammable substances



Safe storage of flammable or toxic gases



Underbench cabinet for safe storage of flammable and corrosive substances



Underbench cabinet with refrigerated unit and recirculating air filter system for the unrestricted storage of flammable hazardous materials in working areas



Underbench cabinet for safe collection of liquid waste



Secure storing and charging solution for lithium-ion batteries



Shock-resistant cabinets for earthquake zones



cabinet for safe storage of oils, varnishes, thinners, colors



Steel furniture manufactured according to European regulations and standards. Numerous options for chemical, microbiological, pharmaceutical, food, R&D control laboratories, as well for production and warehouse facilities. Custom-made configurations of furniture, based on special demands in regard with dimensions, quality of work-tops, end-user's analyses and protocols.





### Model Premium Pro plus

Fume hood of premium quality. Upgraded construction, electrics and electronic control solutions. Comfortable and safe work. Various options regarding materials and supply for all applications.

### Model Premium Pro

Fume hood of highest quality. Full steel frame with chamber lining resistant to chemicals, temperature, and moist. Flexible design for out-of-standard dimensions and applications.



### Model Classic Pro

Fume hood for applications with less thermal and moist load, with high chemical resistance. Ideal for sample preparation and small lab apparatus housing.

### Model Economic Pro

Table top fume hood with lightweight frame and chemical resistant chamber lining. Suitable for easy location change, but still with all supply options.



## WHY CHOOSE US?

Offering custom-made options according to end-users' needs, demands and limits of the area (laboratories, warehouses, production plants, pharmacies, hospitals).

Delivery of furniture, fume hoods and safety cabinets manufactured according to standards EN 16121, EN 14470-1, EN 14470-2, EN 14175.

Highly professional installation of fume hoods and safety cabinets, following local regulations, HSE guidance, offering options of connecting the cabinets to an existing exhaust system or using a recirculating air filter add-on.

After-sales maintenance, according to good-laboratory practice, end-user's protocols and strict recommendations of the manufacturer.

Overall satisfaction of employees in working area and managers, due to fully engaged solutions which provide maximum comfort, top-protection and compliance with EU and local legislation.



## SAFETY MATTERS